

UBE2J2 抗原(重组蛋白)

中文名称: UBE2J2 抗原(重组蛋白)

英文名称: UBE2J2 Antigen (Recombinant Protein)

别 名: ubiquitin-conjugating enzyme E2, J2; NCUBE2; NCUBE-2; PRO2121

相关类别: 抗原

储 存: 冷冻(-20℃)

概述

Fusion protein corresponding to a region derived from 110-259 amino acids of human UBE2J2

技术规格

Full name:	ubiquitin-conjugating enzyme E2, J2
Synonyms:	NCUBE2; NCUBE-2; PRO2121
Swissprot:	Q8N2K1
Gene Accession:	BC113430
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The modification of proteins with ubiquitin is an important cell ular mechanism for targeting abnormal or short-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: ubiquitin-activating enzymes, or E1s, ubiquitin-conjugating enzymes, or E2s, and ubiquitin-protein ligases, or E3s. This gene encodes a member of the E2 ubiquitin-conjugating enzyme family. This enzyme is located in the membrane of the endoplasmic reticulum. Multiple alternatively spliced transcript



variants have been found for this gene, but the full-length nat
ure of some variants has not been defined.