

UBE2Q2 抗原(重组蛋白)

- 中文名称: UBE2Q2 抗原(重组蛋白)
- 英文名称: UBE2Q2 Antigen (Recombinant Protein)
- 别 名: ubiquitin conjugating enzyme E2Q family member 2
- 相关类别: 抗原
- 储 存: 冷冻 (-20℃)

概述

Fusion protein corresponding to a region derived from 176-375 amino acids of human UBE2Q2

技术规格

Full name:	ubiquitin conjugating enzyme E2Q family member 2
Swissprot:	Q8WVN8
Gene Accession:	BC017708
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	UBE2Q2 (Ubiquitin-conjugating enzyme E2 Q2) is a 375 amin o acid cytoplasmic protein that is involved in ubiquitin-media ted protein degradation. Ubiquitination is an important mech anism through which three classes of enzymes act in concert to target short-lived or abnormal proteins for destruction. Th e three classes of enzymes involved in ubiquitination are the ubiquitin-activating enzymes (E1s), the ubiquitin-conjugating enzymes (E2s) and the ubiquitin-protein ligases (E3s). One of



several members of the ubiquitin-conjugating enzyme family,
UBE2Q2 is an E2 ubiquitin-conjugating enzyme that acts to c
atalyze the covalent attachment of ubiquitin residues to vario
us proteins. UBE2Q2 expression is detected in invasive epithe
lial tissue and tumor masses including head and neck squam
ous cell carcinomas, suggesting a role for UBE2Q2 in carcino
genesis. Two isoforms of UBE2Q2 exist due to alternative spli
cing events.