

UBE2U 抗原（重组蛋白）

中文名称： UBE2U 抗原（重组蛋白）

英文名称： UBE2U Antigen (Recombinant Protein)

别名： ubiquitin-conjugating enzyme E2U (putative)

相关类别： 抗原

储存： 冷冻（-20℃）

概述

Fusion protein corresponding to N terminal 226 amino acids of human UBE2U

技术规格

Full name:	ubiquitin-conjugating enzyme E2U (putative)
Swissprot:	Q5VVX9
Gene Accession:	BC029895
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	UBE2U (ubiquitin-conjugating enzyme E2U) is a 321 amino acid E2 ubiquitin conjugating enzyme that catalyzes the covalent attachment of ubiquitin to other proteins. Existing as two alternatively spliced isoforms, UBE2U is encoded by a gene located on human chromosome 1, which spans 260 million base pairs, contains over 3,000 genes and comprises nearly 8% of the human genome.

mlbio 梅联生物
Good elisakit producers