

## UBLCP1 抗原（重组蛋白）

中文名称： UBLCP1 抗原（重组蛋白）

英文名称： UBLCP1 Antigen (Recombinant Protein)

别名： ubiquitin like domain containing CTD phosphatase 1; CPUB1

储存： 冷冻（-20℃）

相关类别： 抗原

### 概述

Fusion protein corresponding to a region derived from 1-318 amino acids of human UBLCP1

### 技术规格

<b>Full name:</b>	ubiquitin like domain containing CTD phosphatase 1
<b>Synonyms:</b>	CPUB1
<b>Swissprot:</b>	Q8WVY7
<b>Gene Accession:</b>	BC013425
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	Dephosphorylates 26S nuclear proteasomes, thereby decreasing their proteolytic activity. The dephosphorylation may prevent assembly of the core and regulatory particles (CP and RP) into mature 26S proteasome.