

USP53 抗原（重组蛋白）

中文名称： USP53 抗原（重组蛋白）

英文名称： USP53 Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 1-188 amino acids of human USP53

技术规格

Full name:	ubiquitin specific peptidase 53
Swissprot:	Q70EK8
Gene Accession:	BC017382
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Tight junction-associated protein that is involved in the survival of a uditory hair cells and hearing. Maybe by modulating the barrier properties and mechanical stability of tight junctions (By similarity). Has no peptidase activity (PubMed:14715245).