

## UTP3 抗原（重组蛋白）

中文名称： UTP3 抗原（重组蛋白）

英文名称： UTP3 Antigen (Recombinant Protein)

别名： UTP3 small subunit processome component; CRL1; CRLZ1; SAS10

相关类别： 抗原

储存： 冷冻（-20℃）

### 概述

Fusion protein corresponding to a region derived from 280-479 amino acids of human UTP3

### 技术规格

<b>Full name:</b>	UTP3 small subunit processome component
<b>Synonyms:</b>	CRL1; CRLZ1; SAS10
<b>Swissprot:</b>	Q9NQZ2
<b>Gene Accession:</b>	BC004546
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	Essential for gene silencing: has a role in the structure of silenced chromatin. Plays a role in the developing brain (By similarity).

**mlbio** 梅联生物  
Good elisakit producers