

## VILL 抗原（重组蛋白）

中文名称： VILL 抗原（重组蛋白）

英文名称： VILL Antigen (Recombinant Protein)

别名： villin-like

相关类别： 抗原

储存： 冷冻（-20℃）

### 概述

Fusion protein corresponding to a region derived from 657-856 amino acids of human VILL

### 技术规格

<b>Full name:</b>	villin-like
<b>Swissprot:</b>	O15195
<b>Gene Accession:</b>	BC004300
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	The protein encoded by this gene belongs to the villin/gelsolin family. It contains 6 gelsolin-like repeats and a headpiece domain. It may play a role in actin-bundling.

**mlbio** 梅联生物  
Good elisakit producers