

VPREB3 抗原(重组蛋白)

- 中文名称: VPREB3 抗原(重组蛋白)
- 英文名称: VPREB3 Antigen (Recombinant Protein)
- 别 名: pre-B lymphocyte 3; 8HS20; N27C7-2
- 相关类别: 抗原
- 储 存: 冷冻(-20℃)

概述

Fusion protein corresponding to a region derived from 21-123 amino acids of human VPREB3

技术规格

Full name:	pre-B lymphocyte 3
Synonyms:	8HS20; N27C7-2
Swissprot:	Q9UKI3
Gene Accession:	BC020666
Purity:	>85%, as determined by Coomassie blue stained SDS-PAG E
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene is the human ortholog of the mouse VpreB3 (8HS20) protein, is thought to be in volved in B-cell maturation, and may play a role in assem bly of the pre-B cell receptor (pre-BCR). While the role of this protein in B-cell development has not yet been elucid ated, studies with the chicken ortholog of this protein hav



e found that when overexpressed, this protein localizes to
the endoplasmic reticulum. The mouse ortholog of this pro
tein has been shown to associate with membrane mu heav
y chains early in the course of pre-B cell receptor biosynth
esis. Expression of this gene has been observed in some ly
mphomas.