

VPS54 抗原（重组蛋白）

中文名称： VPS54 抗原（重组蛋白）

英文名称： VPS54 Antigen (Recombinant Protein)

别名： VPS54 subunit of GARP complex; WR; HCC8; SLP-8p; VPS54L; hVps54L; PPP1R164

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to C terminal 200 amino acids of human VPS54

技术规格

Full name:	VPS54 subunit of GARP complex
Synonyms:	WR; HCC8; SLP-8p; VPS54L; hVps54L; PPP1R164
Swissprot:	Q9P1Q0
Gene Accession:	BC030275
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes for a protein that in yeast forms part of a trimeric vacuolar-protein-sorting complex that is required for retrograde transport of proteins from prevacuoles to the late Golgi compartment. As in yeast, mammalian Vps54 proteins contain a coiled-coil region and dileucine motifs. Alternative splicing results in multiple transcript variants encoding different isoforms.