

WDFY2 抗原（重组蛋白）

中文名称： WDFY2 抗原（重组蛋白）

英文名称： WDFY2 Antigen (Recombinant Protein)

别名： WD repeat and FYVE domain containing ; PROF; WDF2; ZFYVE22

相关类别： 抗原

储存： 冷冻（-20℃）

概述

Fusion protein corresponding to a region derived from 201-400 amino acids of human WDFY2

技术规格

Full name:	WD repeat and FYVE domain containing
Synonyms:	PROF; WDF2; ZFYVE22
Swissprot:	Q96P53
Gene Accession:	BC014004
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a protein that contains two WD domains and an FYVE zinc finger region. The function of this gene is unknown. An alternatively spliced transcript variant of this gene may exist.

mlbio 梅联生物
Good elisakit producers