

WDR1 抗原（重组蛋白）

中文名称：WDR1 抗原（重组蛋白）

英文名称：WDR1 Antigen (Recombinant Protein)

别名：WD repeat domain 1; AIP1; NORI-1; HEL-S-52

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 407-606 amino acids of human WDR1

技术规格

Full name:	WD repeat domain 1
Synonyms:	AIP1; NORI-1; HEL-S-52
Swissprot:	O75083
Gene Accession:	BC002489
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a protein containing 9 WD repeats. WD repeats are approximately 30- to 40-amino acid domains containing several conserved residues, mostly including a trp-asp at the C-terminal end. WD domains are involved in protein-protein interactions. The encoded protein may help induce the disassembly of actin filaments. Two transcript variants encoding different isoforms have been found for this gene.