

WDR18 抗原(重组蛋白)

中文名称: WDR18 抗原(重组蛋白)

英文名称: WDR18 Antigen (Recombinant Protein)

别 名: WD repeat domain 18; Ipi3; R32184_1

储 存: 冷冻(-20℃)

相关类别: 抗原

概述

Fusion protein corresponding to a region derived from 233-432 amino acids of human WDR18

技术规格

Full name:	WD repeat domain 18
Synonyms:	Ipi3; R32184_1
Swissprot:	Q9BV38
Gene Accession:	BC001648
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the WD repeat protein family. WD repeats are minimally conserved regions of approximately 40 amino acids typically bracketed by gly-his and trp-asp (GH-WD), which may facilitate formation of heterotrimeric or multip rotein complexes. Members of this family are involved in a vari ety of cellular processes, including cell cycle progression, signal transduction, apoptosis, and gene regulation.