

WSB2 抗原（重组蛋白）

中文名称： WSB2 抗原（重组蛋白）

英文名称： WSB2 Antigen (Recombinant Protein)

别名： WD repeat and SOCS box containing 2; SBA2

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 255-404 amino acids of human WSB2

技术规格

Full name:	WD repeat and SOCS box containing 2
Synonyms:	SBA2
Swissprot:	Q9NYS7
Gene Accession:	BC015887
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the WD-protein subfamily. The encoded protein contains five WD-repeats spanning most of the protein and an SOCS box in the C-terminus. The SOCS box may act as a bridge between specific substrate-binding domains and E3 ubiquitin protein ligases. Three transcript variants encoding different isoforms have been found for this gene.