

## XRN2 抗原（重组蛋白）

中文名称： XRN2 抗原（重组蛋白）

英文名称： XRN2 Antigen (Recombinant Protein)

别名： 5'-3' exoribonuclease 2

相关类别： 抗原

储存： 冷冻（-20℃）

### 概述

Fusion protein corresponding to a region derived from 751-950 amino acids of human XRN2

### 技术规格

<b>Full name:</b>	5'-3' exoribonuclease 2
<b>Swissprot:</b>	Q9H0D6
<b>Gene Accession:</b>	BC006417
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene shares similarity with the mouse Dhml and the yeast dhp1 gene. The yeast gene is involved in homologous recombination and RNA metabolism, such as RNA synthesis and RNA trafficking. Complementation studies show that Dhml has a similar function in mouse as dhp1. The function of the human gene has not yet been determined. Transcript variants encoding different isoforms have been noted for this gene; however, their full-length nature is not known.

**mlbio** 梅联生物  
Good elisakit producers