

## YY1 抗原（重组蛋白）

中文名称： YY1 抗原（重组蛋白）

英文名称： YY1 Antigen (Recombinant Protein)

别名： YY1 transcription factor; DELTA; NF-E1; UCRBP; INO80S; YIN-YANG-1

储存： 冷冻（-20℃）

相关类别： 抗原

### 概述

Fusion protein corresponding to a region derived from 215-414 amino acids of human YY1

### 技术规格

<b>Full name:</b>	YY1 transcription factor
<b>Synonyms:</b>	DELTA; NF-E1; UCRBP; INO80S; YIN-YANG-1
<b>Swissprot:</b>	P25490
<b>Gene Accession:</b>	BC037308
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	YY1 is a ubiquitously distributed transcription factor belonging to the GLI-Kruppel class of zinc finger proteins. The protein is involved in repressing and activating a diverse number of promoters. YY1 may direct histone deacetylases and histone acetyltransferases to a promoter in order to activate or repress th

e promoter, thus implicating histone modification in the function of YY1.