

ZC3H12A 抗原（重组蛋白）

中文名称： ZC3H12A 抗原（重组蛋白）

英文名称： ZC3H12A Antigen (Recombinant Protein)

别名： zinc finger CCCH-type containing 12A; MCPIP; MCPIP1; dJ423B22.1

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 301-458 amino acids of human ZC3H12A

技术规格

Full name:	zinc finger CCCH-type containing 12A
Synonyms:	MCPIP; MCPIP1; dJ423B22.1
Swissprot:	Q5D1E8
Gene Accession:	BC005001
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	ZC3H12A is an MCP1 (CCL2; MIM 158105)-induced protein that acts as a transcriptional activator and causes cell death of cardiomyocytes, possibly via induction of genes associated with apoptosis.