

ZNF41 抗原(重组蛋白)

- 中文名称: ZNF41 抗原(重组蛋白)
- 英文名称: ZNF41 Antigen (Recombinant Protein)
- 别 名: zinc finger protein 41; MRX89
- 相关类别: 抗原
- 储 存: 冷冻 (-20℃)

概述

Fusion protein corresponding to a region derived from 622-821 amino acids of human ZNF41

技术规格

Full name:	zinc finger protein 41
Synonyms:	MRX89
Swissprot:	P51814
Gene Accession:	BC015023
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene product is a likely zinc finger family transcription f actor. It contains KRAB-A and KRAB-B domains that act as tr anscriptional repressors in related proteins, and multiple zinc finger DNA binding motifs and finger linking regions charact eristic of the Kruppel family. This gene is part of a gene clu ster on chromosome Xp11.23. Several alternatively spliced tr anscript variants have been described, however, the full-lengt



h nature of only some of them is known.