

ZNF43 抗原(重组蛋白)

- 中文名称: ZNF43 抗原(重组蛋白)
- 英文名称: ZNF43 Antigen (Recombinant Protein)
- 别 名: zinc finger protein 43; HTF6; KOX27; ZNF39L1
- 相关类别: 抗原
- 储 存: 冷冻(-20℃)

概述

Fusion protein corresponding to a region derived from 610-809 amino acids of human ZNF43

技术规格

Full name:	zinc finger protein 43
Synonyms:	HTF6; KOX27; ZNF39L1
Swissprot:	P17038
Gene Accession:	BC006528
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene belongs to the C2H2-type zinc finger gene family. Th e zinc finger proteins are involved in gene regulation and devel opment, and are quite conserved throughout evolution. Like this gene product, a third of the zinc finger proteins containing C2H 2 fingers also contain the KRAB domain, which has been found to be involved in protein-protein interactions.

