

ZNF202 抗原（重组蛋白）

中文名称： ZNF202 抗原（重组蛋白）

英文名称： ZNF202 Antigen (Recombinant Protein)

别名： zinc finger protein 202; ZSCAN42; ZKSCAN10

相关类别： 抗原

储存： 冷冻（-20℃）

概述

Fusion protein corresponding to a region derived from 449-648 amino acids of human ZNF202

技术规格

Full name:	zinc finger protein 202
Synonyms:	ZSCAN42; ZKSCAN10
Swissprot:	O95125
Gene Accession:	BC013382
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Transcriptional repressor that binds to elements found predominantly in genes that participate in lipid metabolism. Among its targets are structural components of lipoprotein particles (apolipoproteins AIV, CIII, and E), enzymes involved in lipid processing (lipoprotein lipase, lecithin cholesteryl ester transferase), transporters involved in lipid homeostasis (ABCA1, ABCG1), and several genes involved in processes related to energy metabolism and vascular disease.

mlbio 梅联生物
Good elisakit producers