

ZNF300 抗原（重组蛋白）

中文名称： ZNF300 抗原（重组蛋白）

英文名称： ZNF300 Antigen (Recombinant Protein)

别名： zinc finger protein 300

相关类别： 抗原

储存： 冷冻（-20℃）

概述

Fusion protein corresponding to a region derived from 405-604 amino acids of human ZNF300

技术规格

Full name:	zinc finger protein 300
Swissprot:	Q96RE9
Gene Accession:	BC117248
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene is a C2H2-type zinc finger DNA binding protein and likely transcriptional regulator. The function of this protein is not yet known. Three transcript variants encoding different isoforms have been found for this gene.

mlbio 梅联生物
Good elisakit producers