

人 IgG(H+L)

中文名称： 人 IgG (H+L)

英文名称： Human IgG (H+L)

相关类别： IgG

储 存： 冷藏（2-8℃） 避光

技术规格

Background	Immunoglobulin G (IgG) is one of the most abundant proteins in serum with normal levels between 8-17 mg/mL in adult blood. Ig G is important for our defense against microorganisms and the molecules are produced by B lymphocytes as a part of our adaptive immune response. The IgG molecule has two separate functions , to bind to the pathogen that elicited the response and to recruit other cells and molecules to destroy the antigen. The variability of the IgG pool is generated by somatic recombination and the number of specificities in an individual at a given time point is estimated to be 10 ¹¹ variants.
Applications	Optimal tested applications and working dilutions must be determined by the end user.