

山羊 IgG(H+L)

中文名称: 山羊 IgG (H+L)

英文名称: Goat IgG (H+L)

相关类别: IgG

储 存: 冷冻(-20℃) 避光

技术规格

Background	with normal levels between 8-17 mg/mL in adult blood. IgG is importan t for our defense against microorganisms and the molecules are produc
	ed by B lymphocytes as a part of our adaptive immune response. The I
	gG molecule has two separate functions, to bind to the pathogen that
	elicited the response and to recruit other cells and molecules to destroy
	the antigen. The variability of the IgG pool is generated by somatic rec
	ombination and the number of specificities in an individual at a given ti
	me point is estimated to be 1011 variants.
Applications	Optimal tested applications and working dilutions must be determined
	by the end user.