

生物素标记的人 IgG

中文名称: 生物素标记的人 IgG

英文名称: Biotin-conjugated human IgG

相关类别: IgG

储 存: 冷冻(-20℃) 避光

标 记 物: Biotin

技术规格

Background	Immunoglobulin G (IgG),is one of the most abundant prote ins in serum with normal levels between 9.5-12.5 mg/ml in adult blood. IgG is important for our defense against micr oorganisms and the molecules are produced by B lymphoc ytes as a part of our adaptive immune response. The IgG molecule has two separate functions; to bind to the patho gen that elicited the response and to recruit other cells an d molecules to destroy the antigen. The variability of the I gG pool is generated by somatic recombination and the n umber of specificities in an individual at a given time point is estimated to be 1011 variants.
Applications	1/500-5000 for ELISA 1/500-5000 for Western blotting 1/50 0-5000 for Immunohistochemistry

