

兔抗 DTNB 多克隆抗体

中文名称：兔抗 DTNB 多克隆抗体

英文名称：Anti-DTNB rabbit polyclonal antibody

相关类别：一抗

储存：冷冻（-20℃）

抗原：DTNB

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes dystrobrevin beta, a component of the dystrophin-associated protein complex (DPC). The DPC consists of dystrophin and several integral and peripheral membrane proteins, including dystroglycans, sarcoglycans, syntrophins and dystrobrevin alpha and beta. The DPC localizes to the sarcolemma and its disruption is associated with various forms of muscular dystrophy. Dystrobrevin beta is thought to interact with syntrophin and the DP71 short form of dystrophin. Alternatively spliced transcript variants encoding different isoforms have been identified.
Applications:	ELISA, IHC
Name of antibody:	DTNB
Immunogen:	Fusion protein of human DTNB
Full name:	dystrobrevin, beta
SwissProt:	O60941

ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and Human prostate cancer
IHC Recommend dilution:	50-200

