

## 兔抗 RMC1 多克隆抗体

中文名称: 兔抗 RMC1 多克隆抗体

英文名称: Anti-RMC1 rabbit polyclonal antibody

别名: regulator of MON1-CCZ1; MIC1; Mic-1; WDR98; C18orf

抗原: RMC1

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a colon cancer associated protein. [provided by RefSeq, Jan 2013]
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	RMC1
<b>Immunogen:</b>	Fusion protein of human RMC1
<b>Full name:</b>	regulator of MON1-CCZ1
<b>Synonyms:</b>	MIC1; Mic-1; WDR98; C18orf8; HsT2591
<b>SwissProt:</b>	Q96DM3
<b>IHC positive control:</b>	Human colorectal cancer and Human tonsil
<b>IHC Recommend dilution:</b>	25-100



