

兔抗 HGS 多克隆抗体

中文名称:	兔抗 HGS 多克隆抗体
英文名称:	Anti-HGS rabbit polyclonal antibody
别名:	HRS
相关类别:	一抗
储存:	冷冻 (-20℃) 避光
抗原:	HGS
宿主:	Rabbit
反应种属:	Human Mouse Rat
标记物:	Unconjugate
克隆类型:	Unconjugate

技术规格

Background:	The protein encoded by this gene regulates endosomal sorting and plays a critical role in the recycling and degradation of membrane receptors. The encoded protein sorts monoubiquitinated membrane proteins into the multivesicular body, targeting these proteins for lysosome-dependent degradation.
Applications:	WB, IHC, IF
Name of antibody:	HGS
Immunogen:	Synthesized peptide derived from internal of human HGS.
Full name:	hepatocyte growth factor-regulated tyrosine kinase substrate
Synonyms :	HRS
SwissProt:	O14964
IHC positive control:	Human testis tissue

IHC Recommend dilution:	50-100
WB Predicted band size:	86 kDa
WB Positive control:	293 cells and Jurkat cells lysates
WB Recommended dilution :	500-3000
IF positive control:	Huvec cells
IF Recommend dilution:	100-500



