

## 小鼠抗 PHLDA1 单克隆抗体

中文名称:	小鼠抗 PHLDA1 单克隆抗体
英文名称:	Anti-PHLDA1 mouse monoclonal antibody
别名:	pleckstrin homology like domain family A member 1; PHRIP; TDAG51; DT1P1B11
相关类别:	一抗
储存:	冷冻 (-20℃)
宿主:	Mouse
抗原:	PHLDA1
反应种属:	Human
标记物:	Unconjugate
克隆类型:	mouse monoclonal

### 技术规格

<b>Background:</b>	This gene encodes an evolutionarily conserved proline-histidine rich nuclear protein. The encoded protein may play an important role in the anti-apoptotic effects of insulin-like growth factor-1.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	PHLDA1
<b>Immunogen:</b>	Fusion protein of human PHLDA1
<b>Full name:</b>	pleckstrin homology like domain family A member 1
<b>Synonyms:</b>	PHRIP; TDAG51; DT1P1B11
<b>SwissProt:</b>	Q8WV24
<b>WB Predicted band size:</b>	45 KD
<b>WB Positive control:</b>	A375 cell lysates



**WB Recommended dilution:**

500-2000