

## 兔抗 HOXA2 多克隆抗体

中文名称：兔抗 HOXA2 多克隆抗体

英文名称：Anti-HOXA2 rabbit polyclonal antibody

别名：HOX1K; MCOHI

相关类别：一抗

抗原：HOXA2

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：Unconjugate

### 技术规格

<b>Background:</b>	This gene is part of the A cluster on chromosome 7 and encodes a DNA-binding transcription factor which may regulate gene expression, morphogenesis, and differentiation. The encoded protein may be involved in the placement of hindbrain segments in the proper location along the anterior-posterior axis during development.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	HOXA2
<b>Immunogen:</b>	Synthetic peptide of human HOXA2
<b>Full name:</b>	homeobox A2
<b>Synonyms :</b>	HOX1K; MCOHI
<b>SwissProt:</b>	O43364
<b>ELISA Recommend</b>	5000-10000

<b>ed dilution:</b>	
<b>IHC positive control:</b>	Human colorectal cancer
<b>IHC Recommended dilution:</b>	25-100
<b>WB Predicted band size:</b>	41 kDa
<b>WB Positive control:</b>	Human fetal muscle tissue and PC-3 cell
<b>WB Recommended dilution:</b>	200-1000

