

## 兔抗 IL4 多克隆抗体

中文名称： 兔抗 IL4 多克隆抗体

英文名称： Anti-IL4 rabbit polyclonal antibody

别名： BSF1; IL-4; BCGF1; BSF-1; BCGF-1

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： IL4

反应种属： Human

标记物： Unconjugate


克隆类型： Unconjugate

### 技术规格

**Background:**

The protein encoded by this gene is a pleiotropic cytokine produced by activated T cells. This cytokine is a ligand for interleukin 4 receptor. The interleukin 4 receptor also binds to IL13, which may contribute to many overlapping functions of this cytokine and IL13. STAT6, a signal transducer and activator of transcription, has been shown to play a central role in mediating the immune regulatory signal of this cytokine. This gene, IL3, IL5, IL13, and CSF2 form a cytokine gene cluster on chromosome 5q, with this gene particularly close to IL13. This gene, IL13 and IL5 are found to be regulated coordinately by several long-range regulatory elements in an over 120 kilobase range on the chromosome. Two a

	alternatively spliced transcript variants of this gene encoding distinct isoforms have been reported.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	IL4
<b>Immunogen:</b>	Synthesized peptide derived from C-terminal of human IL4.
<b>Full name:</b>	interleukin 4
<b>Synonyms :</b>	BSF1; IL-4; BCGF1; BSF-1; BCGF-1
<b>SwissProt:</b>	P05112
<b>WB Predicted band size:</b>	17 kDa
<b>WB Positive control:</b>	Hela cells lysate
<b>WB Recommended dilution:</b>	500-3000



**mlbio** 梅联生物  
Good elisakit producers