

## 兔抗 IL1RAPL2 多克隆抗体

中文名称：兔抗 IL1RAPL2 多克隆抗体

英文名称：Anti-IL1RAPL2 rabbit polyclonal antibody

别名：IL1R9; IL-1R9; TIGIRR-1; IL1RAPL-2

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：IL1RAPL2

反应种属：Human


标记物：Unconjugate

克隆类型：Unconjugate

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a member of the interleukin 1 receptor family. This protein is similar to the interleukin 1 accessory proteins, and is most closely related to interleukin 1 receptor accessory protein-like 1 (IL1RAPL1). This gene and IL1RAPL1 are located at a region on chromosome X that is associated with X-linked non-syndromic mental retardation.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	IL1RAPL2
<b>Immunogen:</b>	Synthesized peptide derived from internal of human IL1RAPL2
<b>Full name:</b>	interleukin 1 receptor accessory protein-like 2
<b>Synonyms :</b>	IL1R9; IL-1R9; TIGIRR-1; IL1RAPL-2

<b>SwissProt:</b>	Q9NP60
<b>WB Predicted band size:</b>	79 kDa
<b>WB Positive control:</b>	COLO205 cells lysate
<b>WB Recommended dilution:</b>	500-3000



**mlbio** 酶联生物  
Good elisakit producers

