

兔抗 IGF1 多克隆抗体

中文名称：兔抗 IGF1 多克隆抗体

英文名称：Anti-IGF1 rabbit polyclonal antibody

别名：IGFI, IGF-I, IGF1A

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：IGF1

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is similar to insulin in function and structure and is a member of a family of proteins involved in mediating growth and development. The encoded protein is processed from a precursor, bound by a specific receptor, and secreted. Defects in this gene are a cause of insulin-like growth factor I deficiency. Several transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	IGF1
Immunogen:	Synthetic peptide of human IGF1
Full name:	insulin-like growth factor 1 (somatomedin C)

Synonyms :	IGFI, IGF-I, IGF1A
SwissProt:	P05019
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and Human thyroid cancer
IHC Recommend dilution :	25-100



