

兔抗 ILKAP 多克隆抗体

中文名称:	兔抗 ILKAP 多克隆抗体
英文名称:	Anti-ILKAP rabbit polyclonal antibody
别名:	ILK associated serine/threonine phosphatase; PPM1O; ILKAP2; ILKAP3; PP2C-DELTA
相关类别:	一抗
储存:	冷冻 (-20℃)
抗原:	ILKAP
宿主:	Rabbit
反应种属:	Human, Mouse, Rat
标记物:	Unconjugate
克隆类型:	rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a protein serine/threonine phosphatase of the PP2C family. This protein can interact with integrin-linked kinase (ILK/ILK1), a regulator of integrin mediated signaling, and regulate the kinase activity of ILK. Through the interaction with ILK, this protein may selectively affect the signaling process of ILK-mediated glycogen synthase kinase 3 beta (GSK3beta), and thus participate in Wnt signaling pathway.
Applications:	ELISA, WB, IHC
Name of antibody:	ILKAP
Immunogen:	Fusion protein of human ILKAP
Full name:	ILK associated serine/threonine phosphatase

Synonyms:	PPM1O; ILKAP2; ILKAP3; PP2C-DELTA
SwissProt:	Q9H0C8
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human ovarian cancer
IHC Recommend dilution:	50-300
WB Predicted band size:	43 kDa
WB Positive control:	Hela and 293T cell lysates
WB Recommended dilution:	500-2000



