

兔抗 IQGAP1 多克隆抗体

中文名称: 兔抗 IQGAP1 多克隆抗体

英文名称: Anti-IQGAP1 rabbit polyclonal antibody

级 别: SAR1; p195; HUMORFA01

抗 原: IQGAP1

储 存: 冷冻 (-20℃)

宿 主: Rabbit

反应种属: Human, Mouse

相关类别: 一抗

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

| | |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Background: | This gene encodes a member of the IQGAP family. The protein contains four IQ domains, one calponin homology domain, one Ras-GAP domain and one WW domain. It interacts with components of the cytoskeleton, with cell adhesion molecules, and with several signaling molecules to regulate cell morphology and motility. Expression of the protein is upregulated by gene amplification in two gastric cancer cell lines. |
| Applications: | ELISA, WB, IHC |
| Name of antibody: | IQGAP1 |
| Immunogen: | Synthetic peptide of human IQGAP1 |
| Full name: | IQ motif containing GTPase activating protein 1 |
| Synonyms: | SAR1; p195; HUMORFA01 |

| | |
|------------------------------------|-----------------------|
| SwissProt: | P46940 |
| ELISA Recommended dilution: | 5000-10000 |
| IHC positive control: | Human cervical cancer |
| IHC Recommend dilution: | 25-100 |
| WB Predicted band size: | 189 kDa |
| WB Positive control: | HUVEC cell lysate |
| WB Recommended dilution: | 500-2000 |



