

## 兔抗 HLA-H 多克隆抗体

中文名称: 兔抗 HLA-H 多克隆抗体

英文名称: Anti-HLA-H rabbit polyclonal antibody

别名: HLAHP

抗原: HLA-H

储存: 冷冻 (-20℃) 避光

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

### 技术规格

<b>Background:</b>	This major histocompatibility complex gene represents a transcribed pseudogene, possibly derived from HLA-A. This gene displays extensive variation.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	HLA-H
<b>Immunogen:</b>	Synthesized peptide derived from internal of human HLAH.
<b>Full name:</b>	major histocompatibility complex, class I, H (pseudogene)
<b>Synonyms :</b>	HLAHP
<b>SwissProt:</b>	P01893
<b>WB Predicted band size:</b>	41 kDa
<b>WB Positive control:</b>	LOVO cells lysate
<b>WB Recommended dilution:</b>	500-3000

