

兔抗 HMBOX1 多克隆抗体

中文名称: 兔抗 HMBOX1 多克隆抗体

英文名称: Anti-HMBOX1 rabbit polyclonal antibody

别名: homeobox containing 1; HOT1; PBHNF; HNF1LA

储存: 冷冻 (-20℃)

抗原: HMBOX1 宿主 Rabbit

反应种属: Human, Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Transcription factor. Isoform 1 acts as a transcriptional repressor. Isoform 4 has very low activity as a transcriptional repressor.
Applications:	ELISA, IHC
Name of antibody:	HMBOX1
Immunogen:	Synthetic peptide of human HMBOX1
Full name:	homeobox containing 1
Synonyms:	HOT1; PBHNF; HNF1LA
SwissProt:	Q6NT76
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil and Human cervical cancer
IHC Recommend dilution:	30-150



