

## 兔抗 IRX3 多克隆抗体

中文名称: 兔抗 IRX3 多克隆抗体

英文名称: Anti-IRX3 rabbit polyclonal antibody

别名: iroquois homeobox 3; IRX-1; IRXB1

抗原: IRX3

储存: 冷冻 (-20°C)

宿主: Rabbit

相关类别: 一抗

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	IRX3 is a member of the Iroquois homeobox gene family and plays a role in an early step of neural development. Members of this family appear to play multiple roles during pattern formation of vertebrate embryos.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	IRX3
<b>Immunogen:</b>	Synthetic peptide of human IRX3
<b>Full name:</b>	iroquois homeobox 3
<b>Synonyms:</b>	IRX-1; IRXB1
<b>SwissProt:</b>	P78415
<b>ELISA Recommended diluti</b>	5000-10000

<b>on:</b>	
<b>IHC positive control:</b>	Human esophagus cancer
<b>IHC Recommend dilution:</b>	20-100

