

兔抗 KLF2 多克隆抗体

中文名称：兔抗 KLF2 多克隆抗体

英文名称：Anti-KLF2 rabbit polyclonal antibody

别名：Kruppel like factor 2; LKLF

抗原：KLF2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

| | |
|--------------------|--|
| Background: | The krüppel-type zinc finger transcription factors comprise a conserved family of DNA binding proteins that are important in developmental regulation. The krüppel zinc finger transcription factor was initially identified in <i>Drosophila</i> as a segmentation gene. Krüppel-like factors that have been characterized in mammals include erythroid krüppel-like transcription factor (EKLF), lung krüppel-like transcription factor (LKLF) and gut krüppel-like transcription factor (GKLF). EKLF is expressed primarily in erythroid tissues, and LKLF expression is limited to the lung. GKLF is found predominantly in gut and has been shown to be expressed during growth arrest. In a developi |
|--------------------|--|

| | |
|------------------------------------|--|
| | ng mouse embryo, LKLF is necessary for normal tunica media formation and blood vessel stabilization. LKLF is also sufficient to program quiescence in T cells by negatively regulating the c-Myc-dependent pathway. The gene for human LKLF maps to chromosome 19p13.11. |
| Applications: | ELISA, WB |
| Name of antibody: | KLF2 |
| Immunogen: | Synthetic peptide of human KLF2 |
| Full name: | Kruppel like factor 2 |
| Synonyms: | LKLF |
| SwissProt: | Q9Y5W3 |
| ELISA Recommended dilution: | 5000-10000 |
| WB Predicted band size: | 37 kDa |
| WB Positive control: | A549 cell lysate |
| WB Recommended dilution: | 200-1000 |

