

## 兔抗 KLHL1 多克隆抗体

中文名称: 兔抗 KLHL1 多克隆抗体

英文名称: Anti-KLHL1 rabbit polyclonal antibody

别名: MRP2

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: KLHL1

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>IHC Recommend dilution:</b>	50-200
<b>IHC positive control:</b>	Human thyroid cancer and human brain
<b>SwissProt:</b>	Q9NR64
<b>Synonyms:</b>	MRP2
<b>Full name:</b>	kelch-like family member 1
<b>Immunogen:</b>	Fusion protein of human KLHL1
<b>Name of antibody:</b>	KLHL1
<b>Applications:</b>	ELISA, IHC
<b>Background:</b>	The KLHL1 protein belongs to a family of actin-organizing proteins related to Drosophila Kelch (Nemes et al., 2000 [PubMed 10888605]).



