

## 兔抗 LIN28B 多克隆抗体

中文名称: 兔抗 LIN28B 多克隆抗体

英文名称: Anti-LIN28B rabbit polyclonal antibody

别名: CSDD2; FLJ16517

相关类别: 一抗

储存: 冷冻 (-20℃) 避光

宿主: Rabbit

抗原: LIN28B

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>WB Positive control:</b>	Hela and 293 cells
<b>WB Recommended dilution:</b>	500-1000
<b>WB Predicted band size:</b>	27 kDa
<b>SwissProt:</b>	Q6ZN17
<b>Synonyms:</b>	CSDD2
<b>Full name:</b>	lin-28 homolog B
<b>Immunogen:</b>	Synthetic peptide of human LIN28B
<b>Name of antibody:</b>	LIN28B
<b>Applications:</b>	WB
<b>Background:</b>	The protein encoded by this gene belongs to the lin-28 family, which is characterized by the presence of a cold-shock domain.

in and a pair of CCHC zinc finger domains. This gene is highly expressed in testis, fetal liver, placenta, and in primary human tumors and cancer cell lines. It is negatively regulated by microRNAs that target sites in the 3' UTR, and overexpression of this gene in primary tumors is linked to the repression of let-7 family of microRNAs and derepression of let-7 targets, which facilitates cellular transformation.

