

兔抗 LIPA 多克隆抗体

中文名称：兔抗 LIPA 多克隆抗体

英文名称：Anti-LIPA rabbit polyclonal antibody

别名：lipase A, lysosomal acid type; LAL; CESD

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：LIPA

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes lipase A, the lysosomal acid lipase (also known as cholesterol ester hydrolase). This enzyme functions in the lysosome to catalyze the hydrolysis of cholesteryl esters and triglycerides. Mutations in this gene can result in Wolman disease and cholesteryl ester storage disease. Alternatively spliced transcript variants have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	LIPA
Immunogen:	Fusion protein of human LIPA
Full name:	lipase A, lysosomal acid type
Synonyms:	LAL; CESD
SwissProt:	P38571

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	50-300

