

## 兔抗 LSAMP 多克隆抗体

中文名称: 兔抗 LSAMP 多克隆抗体

英文名称: Anti-LSAMP rabbit polyclonal antibody

别名: LAMP; IGLON3

抗原: LSAMP

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a neuronal surface glycoprotein found in cortical and subcortical regions of the limbic system. During development of the limbic system, this encoded protein is found on the surface of axonal membranes and growth cones, where it acts as a selective homophilic adhesion molecule, and guides the development of specific patterns of neuronal connections.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	LSAMP
<b>Immunogen:</b>	Fusion protein of human LSAMP
<b>Full name:</b>	limbic system-associated membrane protein

<b>Synonyms :</b>	LAMP; IGLON3
<b>SwissProt:</b>	Q13449
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human brain
<b>IHC Recommend dilution:</b>	25-100



