

兔抗 GAB1 (Phospho-Tyr627)多克隆抗体

中文名称:	兔抗 GAB1 (Phospho-Tyr627)多克隆抗体
英文名称:	Anti-GAB1 (Phospho-Tyr627) rabbit polyclonal antibody
相关类别:	一抗
储存:	冷冻 (-20℃) 避光
宿主:	Rabbit
抗原:	GAB1 (Phospho-Tyr627)
反应种属:	Human, Mouse
标记物:	Unconjugate
克隆类型:	rabbit polyclonal

技术规格

Background:	The protein encoded by Gab1 is a member of the IRS1-like multisubstrate docking protein family. It is an important mediator of branching tubulogenesis and plays a central role in cellular growth response, transformation and apoptosis. Two transcript variants encoding different isoforms have been found for this gene.
Applications:	WB, IF
Name of antibody:	GAB1 (Phospho-Tyr627)
Immunogen:	Synthetic peptide of human GAB1 (Phospho-Tyr627)
Full name:	Growth factor receptor bound protein 2-associated binding protein 1 (Phospho-Tyr627)
SwissProt:	Q13480
WB Predicted band size:	77 kDa
WB Positive control:	293 cells treated with EGF
WB Recommended dilution:	500-1000

IF Positive control:	Hela cells
IF Recommended dilution	100-200



