

兔抗 GAB1 多克隆抗体

中文名称: 兔抗 GAB1 多克隆抗体

英文名称: Anti-GAB1 rabbit polyclonal antibody

相关类别: 一抗

抗 原: GAB1

储 存: 冷冻 (-20℃)

宿 主: Rabbit

反应种属: Human, Mouse

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a member of the IRS1-like multisubstrate docking protein family. It is an important mediator of branching tubulogenesis and plays a central role in cellular growth response, transformation and apoptosis. Two transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	GAB1
Immunogen:	Fusion protein of human GAB1
Full name:	GRB2-associated binding protein 1
SwissProt:	Q13480
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human thyroid cancer and human brain
IHC Recommend dilution:	25-100

