

## 兔抗 GALK2 多克隆抗体

中文名称: 兔抗 GALK2 多克隆抗体

英文名称: Anti-GALK2 rabbit polyclonal antibody

别名: galactokinase 2; GK2

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: GALK2

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a highly efficient N-acetylgalactosamine (GalNAc) kinase, which has galactokinase activity when galactose is present at high concentrations. The encoded protein is a member of the GHMP kinase family. Alternative splicing results in multiple transcript variants.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	GALK2
<b>Immunogen:</b>	Fusion protein of human GALK2
<b>Full name:</b>	galactokinase 2
<b>Synonyms:</b>	GK2
<b>SwissProt:</b>	Q01415
<b>ELISA Recommended dilution:</b>	5000-10000

<b>IHC positive control:</b>	Human tonsil and Human cervical cancer
<b>IHC Recommend dilution:</b>	50-300
<b>WB Predicted band size:</b>	50 kDa
<b>WB Positive control:</b>	293T cell lysate
<b>WB Recommended dilution:</b>	1000-5000



