

兔抗 GALR3 多克隆抗体

中文名称：兔抗 GALR3 多克隆抗体

英文名称： Anti-GALR3 rabbit polyclonal antibody

别名： galanin receptor 3

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： GALR3

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The neuropeptide galanin modulates a variety of physiologic processes including cognition/memory, sensory/pain processing, hormone secretion, and feeding behavior. The human galanin receptors are G protein-coupled receptors that functionally couple to their intracellular effector through distinct signaling pathways. GALR3 is found in many tissues and may be expressed as 1.4-, 2.4-, and 5-kb transcripts
Applications:	ELISA, WB, IHC
Name of antibody:	GALR3
Immunogen:	Synthetic peptide of human GALR3
Full name:	galanin receptor 3

SwissProt:	O60755
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human brain and Human esophagus cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	40 kDa
WB Positive control:	Raji cell lysate
WB Recommended dilution:	500-2000



