

## 兔抗 GAST 多克隆抗体

中文名称：兔抗 GAST 多克隆抗体

英文名称： Anti-GAST rabbit polyclonal antibody

别名： GAS

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： GAST

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

|                          |   |
|--------------------------|---|
| <b>Background:</b>       | Gastrin is a hormone whose main function is to stimulate secretion of hydrochloric acid by the gastric mucosa, which results in gastrin formation inhibition. This hormone also acts as a mitogenic factor for gastrointestinal epithelial cells. Gastrin has two biologically active peptide forms, G34 and G17. |
| <b>Applications:</b>     | ELISA, IHC  |
| <b>Name of antibody:</b> | GAST  |
| <b>Immunogen:</b>        | Synthetic peptide of human GAST   |
| <b>Full name:</b>        | gastrin   |
| <b>Synonyms:</b>         | GAS   |
| <b>SwissProt:</b>        | P01350  |

|                                    |  |
|------------------------------------|--|
| <b>ELISA Recommended dilution:</b> | 5000-10000                             |
| <b>IHC positive control:</b>       | Human esophagus cancer and human brain |
| <b>IHC Recommend dilution:</b>     | 20-100                                 |



