

兔抗 GCGR 多克隆抗体

中文名称：兔抗 GCGR 多克隆抗体

英文名称：Anti-GCGR rabbit polyclonal antibody

别名：GGR

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：GCGR

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a glucagon receptor that is important in controlling blood glucose levels. Defects in this gene are a cause of non-insulin-dependent diabetes mellitus (NIDDM).
Applications:	ELISA, WB, IHC
Name of antibody:	GCGR

Immunogen:	Synthetic peptide of human GCGR
Full name:	glucagon receptor
Synonyms :	GGR
SwissProt:	P47871
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	54 kDa
WB Positive control:	293T cells
WB Recommended dilution:	200-1000



