

兔抗 TRUB1 多克隆抗体

中文名称: 兔抗 TRUB1 多克隆抗体

英文名称: Anti-TRUB1 rabbit polyclonal antibody

别名: TruB pseudouridine synthase family member 1; PUS4

储存: 冷冻 (-20℃)

抗原: TRUB1

宿主: Rabbit

相关类别: 一抗

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Pseudouridine is an abundant component of rRNAs and tRNAs and is enzymatically generated by isomerization of uridine by pseudouridine synthase (Zucchini et al., 2003 [PubMed 12736709]).
Applications:	ELISA, IHC
Name of antibody:	TRUB1
Immunogen:	Fusion protein of human TRUB1
Full name:	TruB pseudouridine synthase family member 1
Synonyms:	PUS4
SwissProt:	Q8WWH5
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil

IHC Recommend dilution:	50-300
-------------------------	--------

