

兔抗 TUFT1 多克隆抗体

中文名称：兔抗 TUFT1 多克隆抗体

英文名称：Anti-TUFT1 rabbit polyclonal antibody

别名：tuftelin 1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：TUFT1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Tuftelin is an acidic protein that is thought to play a role in dental enamel mineralization and is implicated in caries susceptibility. It is also thought to be involved with adaptation to hypoxia, mesenchymal stem cell function, and neurotrophin nerve growth factor mediated neuronal differentiation.
Applications:	ELISA, IHC
Name of antibody:	TUFT1
Immunogen:	Fusion protein of human TUFT1
Full name:	tuftelin 1
SwissProt:	Q9NNX1
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human cervical cancer and Human ovarian cancer
IHC Recommend dilution:	100-300



