

## 兔抗 U2AF2 多克隆抗体

中文名称: 兔抗 U2AF2 多克隆抗体

英文名称: Anti-U2AF2 rabbit polyclonal antibody

别名: U2AF65

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: U2AF2

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

**Background:**

U2 auxiliary factor (U2AF), comprised of a large and a small subunit, is a non-snRNP protein required for the binding of U2 snRNP to the pre-mRNA branch site. This gene encodes the U2AF large subunit which contains a sequence-specific RNA-binding region with 3 RNA recognition motifs and an Arg/Ser-rich domain necessary for splicing. The large subunit binds to the polypyrimidine tract of introns early during spliceosome assembly. Multiple transcript variants have been detected for this gene, but the full-length nature of only two have been detected.

	rmind to date.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	U2AF2
<b>Immunogen:</b>	Fusion protein of human U2AF2
<b>Full name:</b>	U2 small nuclear RNA auxiliary factor 2
<b>Synonyms:</b>	U2AF65
<b>SwissProt:</b>	P26368
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human esophagus cancer
<b>IHC Recommend dilution:</b>	25-100





